



# Introducing ARveo

Digital Augmented Reality Microscope



## ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience



# Augmented Reality Imaging

## GLOW800

### ARveo Platform

#### AR Imaging

GLOW AR Technology

**GLOW800**

Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience

### Vascular fluorescence augmented with GLOW800

Observe cerebral anatomy in natural color, augmented by real-time vascular flow, with full depth perception. GLOW800 AR fluorescence takes the high contrast of ICG and combines it with white light.

Touch to see this technology at work

**GLOW800**  
AR Fluorescence

White Light

NIR Fluorescence



Real-time GLOW800 visualization. Natural color anatomy and high-contrast, colored blood flow with full depth perception.

- > No more interrupting surgery to switch between the natural microscope image and a flat black and white NIR video, no more mental gymnastics to recall and reconcile the different views
- > Crisp delineation helps you limit potential compromise or obstruction of surrounding perforators and small vessels.
- > Full depth perception and no dark peripheries supports clear spatial orientation to aid manipulation of vessels.

**“GLOW800 worked straight away. Suddenly we had the blood vessels lighting up but we could still see the brain structures around them. Now we can keep oriented in the surgical environment.”**

**Professor Raphael Guzman MD**

GLOW800 cannot be viewed in 3D

# Augmented Reality Imaging

## GLOW800

### ARveo Platform

#### AR Imaging

GLOW AR Technology

**GLOW800**

Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience

### Vascular fluorescence augmented with GLOW800


Observe cerebral anatomy in natural color, augmented by real-time vascular flow, with full depth perception. GLOW800 AR fluorescence takes the high contrast of ICG and combines it with white light.

Touch to see this technology at work

**GLOW800**  
AR Fluorescence

**White Light**

**NIR Fluorescence**



- > No more interrupting surgery to switch between the natural microscope image and a flat black and white NIR video, no more mental gymnastics to recall and reconcile the different views
- > Crisp delineation helps you limit potential compromise or obstruction of surrounding perforators and small vessels.
- > Full depth perception and no dark peripheries supports clear spatial orientation to aid manipulation of vessels.

**“GLOW800 worked straight away. Suddenly we had the blood vessels lighting up but we could still see the brain structures around them. Now we can keep oriented in the surgical environment.”**

**Professor Raphael Guzman MD**

GLOW800 cannot be viewed in 3D

# Augmented Reality Imaging

## GLOW800

### ARveo Platform

#### AR Imaging

GLOW AR Technology

#### GLOW800

Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience

### Vascular fluorescence augmented with GLOW800

Observe cerebral anatomy in natural color, augmented by real-time vascular flow, with full depth perception. GLOW800 AR fluorescence takes the high contrast of ICG and combines it with white light.

Touch to see this technology at work

**GLOW800**  
AR Fluorescence

**White Light**

**NIR Fluorescence**

Delay with NIR fluorescence. Flat image lacks anatomical detail and has dark peripheries.



- > No more interrupting surgery to switch between the natural microscope image and a flat black and white NIR video, no more mental gymnastics to recall and reconcile the different views
- > Crisp delineation helps you limit potential compromise or obstruction of surrounding perforators and small vessels.
- > Full depth perception and no dark peripheries supports clear spatial orientation to aid manipulation of vessels.

**“GLOW800 worked straight away. Suddenly we had the blood vessels lighting up but we could still see the brain structures around them. Now we can keep oriented in the surgical environment.”**

**Professor Raphael Guzman MD**

GLOW800 cannot be viewed in 3D

# Augmented Reality Imaging

## A Decade of Pioneering Fluorescence

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

#### Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience


Touch to see the options

Fluorescence Options

FL400

FL800

FL560



#### FL400 oncological fluorescence

- > The fluorescence module FL400 is used during open neurosurgery in conjunction with the active substance 5 aminolevulinic acid (5-ALA). It supports resection by allowing differentiation of tumor tissue from healthy brain tissue.

Please check with your local Leica Microsystems representative for product registration status in your country.

# Augmented Reality Imaging

## A Decade of Pioneering Fluorescence

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

#### Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience


Touch to see the options

Fluorescence  
Options

FL400

FL800

FL560



#### FL800 intraoperative video-angiography module

- > When used in conjunction with ICG fluorescent dye, this feature allows neurosurgeons to clearly visualize blood flow in real-time.

VIDEO FL400

VIDEO FL560

# Augmented Reality Imaging

## A Decade of Pioneering Fluorescence

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

#### Pioneering Fluorescence

IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience


Touch to see the options

Fluorescence  
Options

FL400

FL800

FL560



#### FL560

- > The FL560 fluorescence module enables simultaneous, real-time observation of both non-fluorescent tissue and fluorescent areas, with clear differentiation and contrast.

Please check with your local Leica Microsystems representative for product registration status in your country.



# Augmented Reality Imaging

## IGS and microscope navigation

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

Pioneering Fluorescence

#### IGS with Brainlab

Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience

View images from IGS systems to further augment your surgical insight and enjoy robotic control for seamless workflow. Inject data into the oculars with CaptiView or display on a large 4K 3D screen.




Image shows Cranial Navigation from Brainlab. Brainlab is a trademark of Brainlab AG.

#### Cranial and Microscope Navigation from Brainlab

- > Stay focused on your patient with an all-in-one visualization tool
- > Augmented Reality view within the microscope
- > Guiding tools for improved decision making
- > Anatomical and functional information brought into the surgical field
- > Shift correction with reference to anatomical structures
- > Tool tracking through IGS robotic function

GLOW800 cannot be viewed in 3D

Leica Microsystems | ARveo Customer Presentation | MC-0002465-14.04.2021

Leica Microsystems GmbH. Registered Office: Wetzlar. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.

# Augmented Reality Imaging

## CaptiView Image Injection

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

Pioneering Fluorescence

IGS with Brainlab

#### Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library


#### Case Studies

#### The Leica Demo Experience

VIDEO CaptiView Image Injection

Neurosurgical procedures require your complete focus, so we've put you at the core of our design.

**With CaptiView image injection there's no need to look away from the surgical site to the monitor during surgery.**



- > View GLOW AR fluorescence, IGS data from leading manufacturers, endoscope images and microscope information thanks to the OpenArchitecture design
- > Rely on full-HD 1080p resolution and 500:1 contrast
- > Overlay data on the live surgical image or view as non-correlated in left, right or both eyepieces

GLOW800 cannot be viewed in 3D

# Augmented Reality Imaging

## CaptiView Image Injection

### ARveo Platform

#### AR Imaging

GLOW AR Technology

GLOW800

Pioneering Fluorescence

IGS with Brainlab

#### Image Injection

#### 3D in the OR

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library


#### Case Studies

#### The Leica Demo Experience

VIDEO CaptiView Image Injection

Neurosurgical procedures require your complete focus, so we've put you at the core of our design.

**With CaptiView image injection there's no need to look away from the surgical site to the monitor during surgery.**



- > View GLOW AR fluorescence, IGS data from leading manufacturers, endoscope images and microscope information thanks to the OpenArchitecture design
- > Rely on full-HD 1080p resolution and 500:1 contrast
- > Overlay data on the live surgical image or view as non-correlated in left, right or both eyepieces

GLOW800 cannot be viewed in 3D

# 3D in the OR

## Heads-up Surgery

ARveo Platform

AR Imaging

3D in the OR

Heads-up Surgery

Teaching

Choose Your Monitor

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics


Safety

Video Library

Case Studies

The Leica Demo  
Experience

Choose to work heads-up and enjoy the full depth perception and high resolution you require without needing to look through the oculars.



Visualize minute anatomical details on a much larger scale and with natural color differentiation on 31-inch or 55-inch 4K 3D monitors

Feel even more immersed in your surgical procedure thanks to the magnified view and uncompromised working position

Achieve a comfortable upright posture by positioning the optics carrier and cart-mounted 3D monitor exactly where you need them

# 3D in the OR Teaching

ARveo Platform

AR Imaging

**3D in the OR**

Heads-up Surgery

**Teaching**

Choose Your Monitor

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience

While you're focused on the task at hand, the whole team can follow in detail in 2D or 3D on one of the large 4K monitors



3D visualization on screen can help improve your team's understanding of complex cases and may even increase knowledge retention

Shared 3D viewing on large screens also enables additional students and staff to follow your every delicate move, enhancing not only education but also workflow

## 3D in the OR

# Choose Your Monitor

### ARveo Platform

#### AR Imaging

#### 3D in the OR

Heads-up Surgery

Teaching

**Choose Your Monitor**

#### Sharing Beyond the OR

#### Groundbreaking Optics

#### Workflow and Ergonomics

#### Safety

#### Video Library

#### Case Studies

#### The Leica Demo Experience



- > 55-inch 4K 3D cart system
- > 31-inch 4K 3D microscope monitor
- > Cart mounted system for positioning flexibility in your OR



# Sharing Beyond the OR

## See the Difference Beyond the OR

### ARveo Platform

AR Imaging

3D in the OR

### Sharing Beyond the OR

Live Streaming

Recording and Documenting

Full Security

Application

Groundbreaking Optics





Workflow and Ergonomics

Safety

Video Library

Case Studies

The Leica Demo Experience



# Sharing Beyond the OR

## Live Streaming

### ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond the OR

**Live Streaming**

Recording and Documenting

Full Security

Application

Groundbreaking Optics

Workflow and Ergonomics

Safety

Video Library

Case Studies

The Leica Demo Experience



Hospital Network and/or Media Server

📶
Live Stream

**Stream live video\* instantly to mobile or desktop devices for flexible viewing and education beyond the OR with the integrated Med X Change HDMD Pro.**

- > Share your skills live with students and peers
- > Remotely observe your residents as they prepare the surgical site, without leaning over their shoulders

GLOW800 cannot be viewed in 3D

\*Live streaming is not certified for supporting remote operations

Leica Microsystems | ARveo Customer Presentation | MC-0002465-14.04.2021

Leica Microsystems GmbH. Registered Office: Wetzlar. All rights reserved, also regarding any disposal, exploitation, reproduction, editing, distribution, as well as in the event of applications for industrial property rights.



# Sharing Beyond the OR

## Recording and Documenting

### ARveo Platform

AR Imaging

3D in the OR

**Sharing Beyond  
the OR**

Live Streaming

**Recording and  
Documenting**

Full Security

Application

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience



**Hospital Network  
and/or Media Server**



**Capture, Save and Share with Ease**

### Ready to capture and save it all

- > With built-in workflow features, you don't lose time setting up a camera or recording system
- > Customize your recording, editing and video replay with the fully compliant and secure HDMD Pro from Med X Change



### Customize your recording, editing and video replay for later teaching and documentation.

- > Record video and still images to USB or your hospital network via cable or Wifi
- > DICOM/PACS integration allows you to document cases and save with patient data
- > Easily edit your recordings for teaching and peer presentations

# Sharing Beyond the OR Application

## ARveo Platform

AR Imaging

3D in the OR

**Sharing Beyond  
the OR**

Live Streaming

Recording and  
Documenting

Full Security

**Application**

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

Case Studies


The Leica Demo  
Experience



**Hospital Network  
and/or Media Server**

## Application

### MEDXSTREAM® wireless iOS app



- > Add patient information
- > Wi-Fi / Auto transfer of HD still images
- > Create PDF surgical reports with image annotation / mark-ups
- > Audio dictation embedded into PDF surgical reports
- > Export to video editing apps
- > Record 1080p video directly to the MedXStream app via USB 3.0 (lightning connector)

# Groundbreaking Optics

## Small Angle Illumination

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

FusionOptics

Small Angle Illumination

Workflow and  
Ergonomics

Safety

Video Library

Case Studies


The Leica Demo  
Experience

Deeper light penetration, increased detail visibility and improved depth perception.

Touch to see the benefits of this technology

With  
Small Angle Illumination


Without  
Small Angle Illumination



**Small Angle Illumination (SAI)** distributes light more evenly, and reduces shadows in the surgical field, providing: increased detail visibility, improved depth perception, and deeper light penetration.

### Optimal Field of Illumination

Autolris automatically adjusts the diaphragm so that only the visible area is illuminated. This prevents the possibility of drying or burning exposed tissue outside the field of view.



The more you know, the more empowered you become to make the right decisions for your patients. Small Angle Illumination (SAI) combined with bright 400 W Xenon illumination allows light to penetrate to the bottom of deep, narrow cavities.

# Groundbreaking Optics

## Small Angle Illumination

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

**Groundbreaking  
Optics**

FusionOptics

**Small Angle Illumination**

Workflow and  
Ergonomics

Safety

Video Library

Case Studies


The Leica Demo  
Experience

**Deeper light penetration, increased detail visibility and improved depth perception.**

Touch to see the benefits of this technology

With  
Small Angle Illumination


Without  
Small Angle Illumination



**Small Angle Illumination (SAI) distributes light more evenly, and reduces shadows in the surgical field, providing: increased detail visibility, improved depth perception, and deeper light penetration.**

### **Optimal Field of Illumination**

Autolris automatically adjusts the diaphragm so that only the visible area is illuminated. This prevents the possibility of drying or burning exposed tissue outside the field of view.



The more you know, the more empowered you become to make the right decisions for your patients. Small Angle Illumination (SAI) combined with bright 400 W Xenon illumination allows light to penetrate to the bottom of deep, narrow cavities.



# Workflow and Ergonomics

## Workflow that simply flows



### ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

### Workflow and Ergonomics

Robust and Maneuverable

Optimized Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience



### Thoughtfully designed for your comfort and efficiency.

Any giant leap in technology is only useful if it's used. The ARveo has been expertly designed so that it easily adapts to your preferred style of working and body frame. With ergonomics and efficiency factored into every design decision, from software to switch, you can experience all the benefits of augmented reality-enhanced surgery, without interrupting workflow.

# Workflow and Ergonomics New Features

## ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

**Workflow and  
Ergonomics**

**Robust and Maneuverable**


Optimized Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience



# Workflow and Ergonomics

## Safety

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

**Safety**

Video Library

Case Studies


The Leica Demo  
Experience

**BrightCare Plus optimizes the light intensity relative to the working distance.**


Touch to see how these technologies work

**BrightCare Plus**  
Light Intensity


**Autolris**  
Light Diameter



Long working distance.



Decreased working distance at same illumination setting (left) creates burn potential in conventional microscopes.



BrightCare Plus automatically adapts light intensity to the working distance, providing safer illumination (up to 60% reduction of intensity).

# Workflow and Ergonomics

## Safety

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond the OR

Groundbreaking Optics

Workflow and Ergonomics

**Safety**

Video Library

Case Studies


The Leica Demo Experience

**Autolris automatically adjusts the diaphragm so that only the visible area is illuminated.**


Touch to see how these technologies work

**BrightCare Plus**  
Light Intensity


**Autolris**  
Light Diameter



At low magnification, the field of illumination (yellow) fills the field of view (green) completely.



Previously, as magnification increased, the field of view became smaller, but the illumination outside the field of view could potentially cause tissue burns (red).



Autolris automatically works with the zoom, decreasing the field of illumination as the field of view decreases. There is no peripheral illumination to cause tissue burns outside the field of view.



# Video Library

## Video Collection

**ARveo Platform**

**AR Imaging**

**3D in the OR**

**Sharing Beyond  
the OR**

**Groundbreaking  
Optics**


**Workflow and  
Ergonomics**

**Safety**

**Video Library**

**Case Studies**


**The Leica Demo  
Experience**



ARveo  
Introduction




GLOW800  
First Impressions



GLOW800  
Comparison




GLOW800  
AV Fistula




GLOW800  
ICA Aneurysm




GLOW800  
MCA Aneurysm




GLOW800  
Pericallosal Aneurysm



CaptiView  
Image Injection



FL400



FL560



Heads-up surgery  
A new era in neurosurgery

GLOW800 cannot be viewed in 3D

# Case Studies

## GLOW800 Augmented Reality

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics


Workflow and  
Ergonomics


Safety

Video Library

**Case Studies**

The Leica Demo  
Experience

From Eye to Insight 





CASE STUDY  
GLOW800 Augmented Reality Fluorescence  
in Aneurysm Treatment

AUTHORS: **Prof. Dr. Feres Chaddad**  
Head of Vascular Neurosurgery  
Federal University of Sao Paulo, Brazil

**Dr. Robert Ibe**  
Clinical Marketing Manager  
Leica Microsystems, Heerbrugg, Switzerland

**GLOW800 Augmented Reality  
Fluorescence in Aneurysm Treatment**

From Eye to Insight 



CASE STUDY  
GLOW800 Augmented Reality Fluorescence in  
AVM (Arteriovenous Malformation) treatment

AUTHORS: **Prof. Dr. Feres Chaddad**  
Head of Vascular Neurosurgery  
Federal University of Sao Paulo, Brazil

**Dr. Robert Ibe**  
Clinical Marketing Manager  
Leica Microsystems, Heerbrugg, Switzerland

**GLOW800 Augmented Reality  
Fluorescence in AVM Treatment**

## Case Studies

# GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics


Workflow and  
Ergonomics


Safety

Video Library

**Case Studies**

The Leica Demo  
Experience

From Eye to Insight 



CASE STUDY

### GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment


AUTHORS:

Prof. Dr. Feres Chaddad Head of Vascular Neurosurgery Federal University of Sao Paulo, Brazil	Dr. Robert Ibe Clinical Marketing Manager Leica Microsystems, Heerbrugg, Switzerland
---	--


CASE STUDY · GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT 2

**Table of Contents**

I. Initial patient presentation and first surgery.....	3
II. Pre-operative assessment before second surgery.....	4-6
III. Intra-operative course and images.....	7-8
IV. Post surgery.....	9-10
V. Impact of GLOW800 Augmented Reality Fluorescence.....	11



UNIVERSIDADE FEDERAL DE SÃO PAULO  
1933



MICROSYSTEMS

All clinical images and videos supplied for this report are courtesy of Prof. Dr. Feres Chaddad, Head of Vascular Neurosurgery at the Federal University of Sao Paulo (UNIFESP), Sao Paulo, Brazil.

Leica Microsystems would like to thank the University of Sao Paulo for the good cooperation for this case study, and especially Prof. Dr. Feres Chaddad for his work on this report.

[www.neurocirurgiaepm.com.br](http://www.neurocirurgiaepm.com.br) [www.leica-microsystems.com](http://www.leica-microsystems.com)

NEXT

GO BACK

# Case Studies

## GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond the OR

Groundbreaking Optics

Workflow and Ergonomics

Safety

Video Library

Case Studies

The Leica Demo Experience

CASE STUDY | GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT

3

### I. INITIAL PATIENT PRESENTATION AND FIRST SURGERY

**Patient:**

A 58-year-old female, with a medical history of well controlled hypertension presented with confusion and complete left hemiparesis.

**Diagnosis:**

- > A ruptured right MCA M1 aneurysm
- > Unruptured left MCA (Middle Cerebral Artery) and PCOM (Posterior Communicating Artery) aneurysms
- > Her only medications were antihypertensives.

**Treatment decision:**

The decision was made to urgently treat the ruptured right MCA aneurysm with microsurgical clipping, in order to protect the patient from rebleeding from the unstable aneurysm.

The unruptured aneurysms were on the contralateral side and needed another craniotomy. Therefore we opted to wait for full recovery from the first surgery and subsequently treat the unruptured aneurysms with microsurgical clipping.

**Outcome first surgery: Treatment of the ruptured right MCA with microsurgical clipping**

The surgery on the ruptured right MCA aneurysm was successful and she was discharged 24 days post op. At the time of her discharge she had recovered the strength in her left leg and arm, but the left facial weakness remained. The cardiologist that evaluated her during hospitalization maintained her previous medications, with good control of her hypertension.




Image courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

PREV

CASE STUDY | GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT


4

### II. PRE-OPERATIVE ASSESSMENT BEFORE SECOND SURGERY

Three months after the first surgery to treat the right-sided ruptured MCA, the patient returned for microsurgical clipping of the unruptured left MCA and PCOM aneurysms. Her physical examination at this time showed her to have a good general appearance, no cardiac abnormalities or signs of peripheral vascular compromise, or other signs of disease. Her neurological examination showed persisting left central facial weakness but no other neurological abnormalities, including of the cranial nerves.

**Pre-operative images**

Fig. 2: Pre-Operative DSA showing the unruptured aneurysms and the clips from the prior surgery.



**Pre-operative images**

Fig. 3: Showing the posterior circulation. Since we had a PCOM aneurysm to treat, we needed to know if the posterior circulation was dependent on the ICA (a persistent fetal circulation).




Image courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

NEXT

GO BACK

## Case Studies

# GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

**Case Studies**

The Leica Demo  
Experience


PREV

CASE STUDY | GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT 5

The following images show the injection of the left ICA. We can see both aneurysms in the 3D reconstruction. Since we had both a PCOM and an MCA aneurysm to treat, we decided on a Transylvian approach, through a Pterional craniotomy.

**Pre-operative images**

Fig. 4:  
3D reconstruction showing the MCA aneurysm.



**Pre-operative images**

Fig. 5:  
3D reconstruction showing the MCA aneurysm.





Image courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

CASE STUDY | GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT 6

**Pre-operative images**

Fig. 6:  
3D reconstruction showing both the MCA and PCOM aneurysms.



**Pre-operative images**

Fig. 7:  
3D reconstruction showing both the MCA and PCOM aneurysms.




Image courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

NEXT

GO BACK

## Case Studies

# GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

Case Studies

The Leica Demo  
Experience

CASE STUDY · GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT

7

### III. INTRA-OPERATIVE COURSE AND IMAGES

We performed a Pterional craniotomy with the withdrawal of the lesser sphenoid wing. With this craniotomy, we exposed the superior and part of the middle temporal gyrus, the Sylvian fissure, and the inferior and part of the middle frontal gyrus. For the approach, we opened Sylvian fissure, separating the frontal lobe from the temporal lobe, where we could find the MCA with its branches (and the aneurysm on the M1 segment). On the proximal M1 MCA, we entered the carotid cistern, where we could find the MCA and the PCOM aneurysms. Prof. Feres Chaddad opened the Sylvian fissure using sharp dissection, maintaining the dissection on the arachnoid space.

#### Pre-Clipping Videos (click to watch online)

##### Identification of MCA & PCOM aneurysms with GLOW800 Augmented Reality Fluorescence from Leica Microsystems


Before clipping, GLOW800 helped to see the flow running through all the vessels and the aneurysms. GLOW800 also helped to distinguish the areas that were thrombosed, on the MCA aneurysm.



Fig. 8: GLOW800: Pre-Clipping. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

After the recognition of the main structures – ICA, ACA (A1), optic nerve, PCOM, choroidal arteries, and MCA (M1 and M2 branches) – the PCOM aneurysm was clipped. It is vital to be sure that both PCOM and the anterior choroidal arteries are preserved.

#### PCOM Aneurysm clipping (click to watch online)



#### Visual confirmation of the clip position (click to watch)



Fig. 9-10: PCOM clipping & visual confirmation of PCOM clipping. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

CASE STUDY · GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT

8

We then performed the occlusion of the MCA aneurysm. When this aneurysm was clipped we had to be sure that the superior and inferior trunks of M1 did not have a kink or a torsion that could lower the flow in these vessels.

#### MCA Aneurysm during clipping (click to watch online)



#### Visual confirmation of the clip position (click to watch)



Fig. 11-12: MCA clipping & visual confirmation of MCA clipping. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

#### Post-Clipping Videos (click to watch online)

##### Visual confirmation of clip position with GLOW800 Augmented Reality Fluorescence from Leica Microsystems

Post clipping, GLOW800 was fundamental to show the patency of the PCOM and the anterior choroidal arteries, as well as the superior and inferior trunks of MCA. GLOW800 is an instant tool to check the flow distal to the clip, which allows one to continue to work during its acquisition, even positioning the clip.




Fig. 13: GLOW800: Post-clipping video. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNIFESP, Sao Paulo, Brazil.

PREV

NEXT

GO BACK

## Case Studies

# GLOW800 Augmented Reality Fluorescence in Aneurysm Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics

Safety

Video Library

**Case Studies**

The Leica Demo  
Experience

PREV

CASE STUDY · GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT 9

**IV. POST SURGERY**  
**CT showing the clip positions**




Fig. 14: CT showing clip positions. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNFESP, Sao Paulo, Brazil.

**CT Angiogram showing the clip positions**




Fig. 15: CT Angiograms showing clip positions. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNFESP, Sao Paulo, Brazil.

CASE STUDY · GLOW800 AUGMENTED REALITY IN ANEURYSM TREATMENT 10

**No flow into the clipped aneurysms confirmed on the 3D reconstruction**




Fig. 16: 3D-reconstruction image. Courtesy of Prof. Dr. Chaddad, Head of Vascular Neurosurgery, UNFESP, Sao Paulo, Brazil.

**Post Surgery Examination and Discharge**

- > GCS 15, pupils equal, round, and reactive to light
- > Left central facial weakness
- > Normal motor function with muscle strength 5/5 bilaterally
- > She was discharged 4 days post surgery

NEXT

GO BACK

# Case Studies

## GLOW800 Augmented Reality Fluorescence in AVM Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond the OR

Groundbreaking Optics


Workflow and Ergonomics

Safety

Video Library

**Case Studies**

The Leica Demo Experience

From Eye to Insight 



**CASE STUDY**


**GLOW800 Augmented Reality Fluorescence in AVM (Arteriovenous Malformation) treatment**


**AUTHORS:**

<b>Prof. Dr. Feres Chaddad</b> Head of Vascular Neurosurgery Federal University of Sao Paulo, Brazil	<b>Dr. Robert Ibe</b> Clinical Marketing Manager Leica Microsystems, Heerbrugg, Switzerland
--	---

**Table of Contents**

1. Initial patient presentation .....	3
2. Pre-operative assessment .....	4-7
3. Intra-operative course and image .....	8-10
4. Post-surgery .....	10-11
5. Impact of GLOW800 Augmented Reality Fluorescence .....	12





All clinical images and videos supplied for this report are courtesy of Prof. Dr. Feres Chaddad, Head of Vascular Neurosurgery at the Federal University of Sao Paulo (UNIFESP), Sao Paulo, Brazil.  
www.neurocirurgiaepm.com.br

Leica Microsystems would like to thank the University of Sao Paulo for the good cooperation for this case study, and especially Prof. Dr. Feres Chaddad for his work on this report.  
www.leica-microsystems.com

NEXT

GO BACK



# Case Studies

## GLOW800 Augmented Reality Fluorescence in AVM Treatment

### ARveo Platform

### AR Imaging

### 3D in the OR

### Sharing Beyond the OR

### Groundbreaking Optics

### Workflow and Ergonomics

### Safety

### Video Library

### Case Studies

### The Leica Demo Experience

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT

3

#### 1. Initial Patient Presentation

##### Patient:

A 24 year old male patient with a four month history of

- > Multiple seizures, two of which were described as tonic-clonic seizures. The others were focal onset seizures, characterized by loss of consciousness, rapidly followed by left hand automatisms, collapses, and return to consciousness after a few minutes
- > Impaired awareness
- > Severe headaches, especially after the seizures

Neurological examination showed slightly decreased attention and memory function, compared to healthy people of similar age and education. No abnormalities of motor, gait, balance, or sensory function were noted.

The patient was treated with carbamazepine (200 mg, 3 x day), which did not improve the seizures. Further assessments and imaging were performed to investigate an underlying cause (see Figures 1-8).

##### Diagnosis:

- > Left temporal pole AVM: Spetzler-Martin Grade II
- > Nidus size 2.5 cm
- > Multiple supply sources: mostly from the middle cerebral and posterior cerebral arteries
- > Multiple drainage points including superior anastomatic vein to the superior sagittal sinus, and tributaries to the inferior petrosal sinus, basal vein, and great cerebral vein

##### Treatment Decision:

At multidisciplinary review, the decision was made to surgically resect the AVM due to the following reasons:

- > Patient had poor response to anticonvulsant treatment
- > Risk of hemorrhage from this malformation remains high if not resected
- > Treatment of choice for Spetzler-Martin Grade II AVMs is microsurgical resection, because this can provide immediate cure with low risk of a recurrence

The aim of the surgery is to achieve control of the seizures and to prevent a hemorrhagic episode that could lead to a neurological sequelae.

PREV

4

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT

#### 2. Pre-operative Assessment

##### Pre-operative Imaging




Figure 1: T1-weighted brain MRI showing the AVM nidus at the left temporal pole with a venous aneurysm (red arrows).



Figure 2: T2-weighted brain MRI showing the AVM nidus located at left temporal pole, with classic flow voids showing blood flow. The lesion induces left hippocampus edema (hyperintensity signal, red arrow).

NEXT

GO BACK

## Case Studies

# GLOW800 Augmented Reality Fluorescence in AVM Treatment

### ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond  
the OR

Groundbreaking  
Optics

Workflow and  
Ergonomics



Safety

Video Library

**Case Studies**

The Leica Demo  
Experience

PREV



NEXT

GO BACK

# Case Studies

## GLOW800 Augmented Reality Fluorescence in AVM Treatment

ARveo Platform

AR Imaging

3D in the OR

Sharing Beyond the OR

Groundbreaking Optics

Workflow and Ergonomics

Safety

Video Library

**Case Studies**

The Leica Demo Experience

PREV

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT 7

**Pre-operative Assessment Images - Steal**




Figure 7: AP projection. It is possible to observe a vascular steal from the contralateral carotid system (red arrow) as well as from the ipsilateral anterior cerebral artery (yellow arrow).

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT 8

**3. Intra-operative Course and Images**

A pretemporal craniotomy was performed to take advantage of the transylvian, lateral subfrontal, subtemporal and temporoportal corridors. With this craniotomy the superior, middle and inferior temporal gyrus, the Sylvian fissure, the inferior and part of the middle frontal gyrus were exposed.

**AVM First Assessment**






Figure 8: After opening the dura the AVM was evident on the temporal lobe and was also assessed using GLOW800.

**GLOW800 Initial Assessment**



Video 1: AVM First Assessment with GLOW800. Surgeon able to manipulate and explore with stereoscopic view while simultaneously assessing blood flow.



Video 2: AVM assessed on deeper dissection - clearly showing filling during arterial phase. During surgery GLOW800 shows the presence of fistulas that need to be disconnected and their relation to the AVM midus, surrounding vessels and brain parenchyma.

NEXT

GO BACK

## Case Studies

# GLOW800 Augmented Reality Fluorescence in AVM Treatment

### ARveo Platform

### AR Imaging

### 3D in the OR

### Sharing Beyond the OR

### Groundbreaking Optics

### Workflow and Ergonomics

### Safety

### Video Library

### Case Studies

### The Leica Demo Experience

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT

9

#### AVM Resection

GLOW800 helped to understand the superficial configuration of the AVM and to confirm the surgical strategy that was decided earlier.


Following the initial assessment of the AVM, Prof. Chaddad opened the Sylvian fissure using sharp dissection in the arachnoid space. He separated the frontal lobe from the temporal lobe in order to gain a better exposure of the AVM and expose the middle cerebral artery, which was the principal feeder of the AVM.

GLOW800 showed the flow through the vessels and the brain


parenchyma in real time, which allowed to continue working while assessing the anatomical structures. The superficial feeders and shunts were then disconnected from the AVM.

Deeper circumferential dissection was then done, coagulating and cutting off feeder vessels in order to progressively decrease the flow through the AVM.

The branches of the posterior cerebral artery were accessed through the subtemporal corridor.



Video 3: Clipping of Fistula under GLOW800. This shows GLOW800 tolerating movement that aids uninterrupted workflow.



Video 4: Delayed filling of the AVM as the fistulas and feeders are excluded. This helped the surgeon better assess treatment as the surgery progressed.



Video 5: Delayed filling of the AVM as more feeders were excluded prior to total exclusion from the normal circulation.




Video 6: AVM excluded from circulation but not yet excised. No flow into the nidus of the AVM.

10

CASE STUDY > GLOW800 AUGMENTED REALITY IN AVM TREATMENT

#### GLOW800 Post AVM Exclusion



Video 7: GLOW800 after the AVM was excised. GLOW800 confirmed complete exclusion from the normal circulation, patency of the vessels proximal to the lesion and absence of any abnormal early drainage veins.

#### 4. Post Surgery




Figure 8: T1-weighted MRI images showing the complete resection of AVM nidus (red arrows).

PREV

NEXT

GO BACK

## ПРЕДСТАВИТЕЛСТВО И СЕРВИЗ ЗА БЪЛГАРИЯ



**РСР ЕООД**

**Централен офис: гр. София, ул. "Виктор Григорович" 3**

**E-mail : [office@rsr-bg.com](mailto:office@rsr-bg.com)**

**Тел: 0895200109**